

Connector™ Luer Lock | REF C45



Connector™ Luer Lock

The Connector is attached to IV lines for closed administration of hazardous drugs. Connector Luer Lock C45 mates with the PhaSeal Injector and provides sealed drug transfer via injection or infusion.

BILL OF MATERIALS

Pos.	Components	Material	Characteristics	Qty
1	Connector	Polypropylene (PP)	White	1
2	Membrane	Thermoplastic elastomer (TPE)	White	1
3	Protective cover	Polyethylene (PE)	Blue	1

PHYSICAL CHARACTERISTICS

REF / Art. no.	External fitting	Dimension length	Dimension width	Connector internal volume	Packaging: User box	Packaging: Transport box
C45	Luer Lock	34 mm (1.34 in.)	17 mm (0.67 in.)	0.1 ml	50 pcs	200 pcs

ASSEMBLY DRAWING

