

Infusion Adapter | REF C100



Infusion Adapter

The Infusion Adapter C100 provides a dry connection to the IV bag when the drug is administered via an IV pump set. The adapter allows nurses the opportunity to select the tubing they need to spike into the bag. The Infusion Adapter has an in-built PhaSeal Connector to facilitate closed addition to IV bags. A non-vented IV set should be used together with the C100. If a vented set is used, make sure that the air inlet is closed throughout the infusion procedure. Latex- and PVC-free.

BILL OF MATERIALS

Pos.	Components	Material	Characteristics	Qty
1	Protective cap	PP	White	1
2	Spike housing	ABS	Natural opaque	1
3	Membrane	Thermoplastic elastomer (TPE)	Natural opaque	1
4	Spike port	Thermoplastic elastomer (TPE)	Natural opaque	1

PHYSICAL CHARACTERISTICS

REF / Art. no.	Length	Width	Height	Internal volume	User box	Transport box
C100	90 mm	24 mm	43 mm	0.14 ml	50 pcs	200 pcs

ASSEMBLY DRAWING

