

Protectors

The only universally compatible closed-system drug transfer device designed for all standard size drug vials.



P28
Neck = 28mm
Chamber = 60ml
Capacity

P50
Neck = 20mm
Chamber = 50ml
Capacity

P21
Neck = 20mm
Chamber = 20ml
Capacity

P14
Neck = 13mm
Chamber = 20ml
Capacity

Carmel Pharma, Inc. 7029 Huntley Road Suite O Columbus, OH 43229 USA
Phone: 614-841-0504 Toll Free: 866-487-9250 Fax: 614-841-0525 E-mail: info@carmelpharmausa.com

See instructions for use for additional information. SW 35471 2010-02 U.S.