

Secondary Set | REF C61



Secondary Set

The Secondary Set C61 is a standard secondary set for attachment to primary line. It has an in-built PhaSeal Connector, providing a dry connection to facilitate closed addition to IV bags. The C61 has a Flow Stop for easy priming.

BILL OF MATERIALS

Pos.	Components	Material	Characteristics	Qty
1	Protective cap	Polyethylene (PE)	Opaque	1
2	Spike	ABS	White	1
3	Connector	ABS	White	1
4	Membrane	Thermoplastic elastomer (TPE)	White	1
5	Tubing	PVC, DEHP free	Clear	1
6	Clamp	Polypropylene (PP)	White	1
7	Rotating collar	ABS	White	1
8	Male Luer Lock	ABS	Clear	1
9	Cap for Luer Lock adapter with hydrophobic filter	Polyethylene (PE)/Acrylic copolymer	Opaque	1

PHYSICAL CHARACTERISTICS

REF / Art. no.	Dimension length	Priming volume	Connector internal volume	Packaging: User box	Packaging: Transport box
C61	49 cm (19,3 in.)	3.0 ml	0.15 ml	30 pcs	120 pcs

ASSEMBLY DRAWING

