

Syringe Tray | REF M15



Syringe Tray


The Syringe Tray is used for carrying filled syringes during transport to the administration area. It is designed to accommodate 1–30 mL syringes with Injector Luer Lock N35. Non-sterile.



Optional: An additional tray can be used as a cover.

ASSEMBLY DRAWING

M15



BILL OF MATERIALS

Position	Components	Material	Characteristics	Quantity
1	Tray	APET	Transparent	1 pcs

PHYSICAL CHARACTERISTICS

REF / Art. no.	Dimension length	Dimension width	Dimension height	Packaging: Transport box
M15	306 mm	70 mm	45 mm	120 pcs

Carmel Pharma AB Aminogatan 30 SE-431 53 Mölndal Sweden
 Tel: +46 31 703 04 00 Fax: +46 31 703 04 04 E-mail: info@carmelpharma.com Web: www.carmelpharma.com

Consult instructions for use: SW 34.871, 2010-01.

